

## Product datasheet for **GA100681**

### Human CD40L (CD40LG) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CD40LG CRISPRa kit - CRISPR gene activation of human CD40 ligand
Symbol:	CD40LG
UniProt ID:	<a href="#">P29965</a>
Synonyms:	CD40L; CD154; gp39; hCD40L; HIGM1; IGM; IMD3; T-BAM; TNFSF5; TRAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)