

Product datasheet for **GA100685**

Human CD53 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | CD53 CRISPRa kit - CRISPR gene activation of human CD53 molecule |
| Symbol: | CD53 |
| UniProt ID: | <u>P19397</u> |
| Synonyms: | MOX44; TSPAN25 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)