

Product datasheet for GA100702

Human LRBA activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRBA CRISPRa kit - CRISPR gene activation of human LPS responsive beige-like anchor protein
Symbol:	LRBA
UniProt ID:	P50851
Synonyms:	BGL; CDC4L; CVID8; LAB300; LBA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)