

Product datasheet for **GA100715**

Human P cadherin (CDH3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CDH3 CRISPRa kit - CRISPR gene activation of human cadherin 3
Symbol:	CDH3
UniProt ID:	P22223
Synonyms:	CDHP; HJMD; PCAD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)