

Product datasheet for **GA100731**

Human CDK3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CDK3 CRISPRa kit - CRISPR gene activation of human cyclin dependent kinase 3
Symbol:	CDK3
UniProt ID:	Q00526
Synonyms:	cyclin-dependent kinase 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)