

Product datasheet for **GA100741**

Human p16INK4A (CDKN2A) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | CDKN2A CRISPRa kit - CRISPR gene activation of human cyclin dependent kinase inhibitor 2A |
| Symbol: | CDKN2A |
| UniProt ID: | P42771 |
| Synonyms: | ARF; CDK4I; CDKN2; CMM2; INK4; INK4A; MLM; MTS-1; MTS1; P14; P14ARF; P16; P16-INK4A; P16INK4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)