

## Product datasheet for **GA100743**

### Human p18 INK4c (CDKN2C) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CDKN2C CRISPRa kit - CRISPR gene activation of human cyclin dependent kinase inhibitor 2C
Symbol:	CDKN2C
UniProt ID:	<a href="#">P42773</a>
Synonyms:	INK4C; p18; p18-INK4C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)