

Product datasheet for **GA100779**

Human CEACAM3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | CEACAM3 CRISPRa kit - CRISPR gene activation of human carcinoembryonic antigen related cell adhesion molecule 3 |
| Symbol: | CEACAM3 |
| UniProt ID: | <u>P40198</u> |
| Synonyms: | CD66D; CEA; CGM1; W264; W282 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)