

## Product datasheet for **GA100793**

### Human CHES1 (FOXN3) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | FOXN3 CRISPRa kit - CRISPR gene activation of human forkhead box N3  |
| Symbol:       | FOXN3  |
| Components:   | GA100793G1, CHES1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA100793G2, CHES1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA100793G3, CHES1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">O00409</a></u>  |
| Synonyms:     | C14orf116; CHES1; PRO1635  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)  |

#### Product images:



[View online »](#)