

## Product datasheet for **GA100808**

### Human Muscarinic Acetylcholine Receptor M3 (CHRM3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CHRM3 CRISPRa kit - CRISPR gene activation of human cholinergic receptor muscarinic 3
Symbol:	CHRM3
UniProt ID:	<a href="#">P20309</a>
Synonyms:	EGBRS; HM3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)