

Product datasheet for **GA100809**

Human Muscarinic Acetylcholine Receptor M4 (CHRM4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CHRM4 CRISPRa kit - CRISPR gene activation of human cholinergic receptor muscarinic 4
Symbol:	CHRM4
UniProt ID:	P08173
Synonyms:	HM4; M4R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)