

Product datasheet for **GA100813**

Human **CHRNA3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CHRNA3 CRISPRa kit - CRISPR gene activation of human cholinergic receptor nicotinic alpha 3 subunit
Symbol:	CHRNA3
UniProt ID:	P32297
Synonyms:	LNCR2; NACHRA3; PAOD2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)