

Product datasheet for **GA100819**

Human **CHRNB3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CHRNB3 CRISPRa kit - CRISPR gene activation of human cholinergic receptor nicotinic beta 3 subunit
Symbol:	CHRNB3
UniProt ID:	Q05901
Synonyms:	acetylcholine receptor, neuronal nicotinic, beta-3 subunit; cholinergic receptor, nicotinic, beta 3; cholinergic receptor, nicotinic, beta polypeptide 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)