


Product datasheet for **GA100839**

Human AP3S1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AP3S1 CRISPRa kit - CRISPR gene activation of human adaptor related protein complex 3 subunit sigma 1
Symbol:	AP3S1
UniProt ID:	Q92572
Synonyms:	adaptor-related protein complex 3, sigma 1 subunit; CLAPS3; CLAPS3, Sigma3A; clathrin-associated/assembly/adaptor protein, small 3 (22kD); Sigma3A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)