

Product datasheet for **GA100843**

Human **CLCN2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CLCN2 CRISPRa kit - CRISPR gene activation of human chloride voltage-gated channel 2
Symbol:	CLCN2
UniProt ID:	P51788
Synonyms:	CIC-2; cIC-2; CLC2; ECA2; ECA3; EGI3; EGI11; EGMA; EJM6; EJM8; LKPAT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)