

Product datasheet for **GA100844**

Human **CLCN3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CLCN3 CRISPRa kit - CRISPR gene activation of human chloride voltage-gated channel 3
Symbol:	CLCN3
UniProt ID:	P51790
Synonyms:	CIC-3; CLC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)