

Product datasheet for **GA100845**

Human **CLCN4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CLCN4 CRISPRa kit - CRISPR gene activation of human chloride voltage-gated channel 4
Symbol:	CLCN4
UniProt ID:	P51793
Synonyms:	CIC-4; CIC-4A; CLC4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)