

Product datasheet for **GA100847**

Human **CLCN6** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CLCN6 CRISPRa kit - CRISPR gene activation of human chloride voltage-gated channel 6
Symbol:	CLCN6
UniProt ID:	P51797
Synonyms:	CLC-6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)