

## Product datasheet for **GA100858**

### Human CLN3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CLN3 CRISPRa kit - CRISPR gene activation of human CLN3 lysosomal/endosomal transmembrane protein, battenin
Symbol:	CLN3
UniProt ID:	<a href="#">Q13286</a>
Synonyms:	BTN1; BTS; JNCL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)