

Product datasheet for **GA100868**

Human CCR3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CCR3 CRISPRa kit - CRISPR gene activation of human C-C motif chemokine receptor 3
Symbol:	CCR3
UniProt ID:	P51677
Synonyms:	CC-CKR-3; CD193; CKR3; CMKBR3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)