

Product datasheet for **GA100869**

Human CCR4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CCR4 CRISPRa kit - CRISPR gene activation of human C-C motif chemokine receptor 4
Symbol:	CCR4
UniProt ID:	P51679
Synonyms:	CC-CKR-4; CD194; ChemR13; CKR4; CMKBR4; HGCN:14099; K5-5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)