

Product datasheet for **GA100870**

Human CCR5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CCR5 CRISPRa kit - CRISPR gene activation of human C-C motif chemokine receptor 5 (gene/pseudogene)
Symbol:	CCR5
UniProt ID:	P51681
Synonyms:	CC-CKR-5; CCCKR5; CCR-5; CD195; CKR-5; CKR5; CMKBR5; IDDM22
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)