

## Product datasheet for **GA100876**

### Human Leukotriene B4 Receptor (LTB4R) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LTB4R CRISPRa kit - CRISPR gene activation of human leukotriene B4 receptor
Symbol:	LTB4R
UniProt ID:	<a href="#">Q15722</a>
Synonyms:	BLT1; BLTR; CMKRL1; GPR16; LTB4R1; LTBR1; P2RY7; P2Y7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)