

Product datasheet for **GA100881**

Human CNG3 (CNGA3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CNGA3 CRISPRa kit - CRISPR gene activation of human cyclic nucleotide gated channel subunit alpha 3
Symbol:	CNGA3
UniProt ID:	Q16281
Synonyms:	ACHM2; CCNC1; CCNCa; CCNCalpha; CNGG3; CNG3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)