

Product datasheet for **GA100883**

Human PLK3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLK3 CRISPRa kit - CRISPR gene activation of human polo like kinase 3
Symbol:	PLK3
UniProt ID:	Q9H4B4
Synonyms:	CNK; FNK; PLK-3; PRK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)