

Product datasheet for **GA100887**

Human CNPase (CNP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CNP CRISPRa kit - CRISPR gene activation of human 2',3'-cyclic nucleotide 3' phosphodiesterase
Symbol:	CNP
UniProt ID:	P09543
Synonyms:	CNP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)