

Product datasheet for **GA100896**

Human Collagen III (COL3A1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	COL3A1 CRISPRa kit - CRISPR gene activation of human collagen type III alpha 1 chain
Symbol:	COL3A1
UniProt ID:	P02461
Synonyms:	EDS4A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)