

Product datasheet for **GA100899**

Human Collagen IV (COL4A3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	COL4A3 CRISPRa kit - CRISPR gene activation of human collagen type IV alpha 3 chain
Symbol:	COL4A3
UniProt ID:	Q01955
Synonyms:	alpha 3 type IV collagen; alpha3 type IV collagen; collagen, type IV, alpha 3 (Goodpasture antigen); collagen IV, alpha-3 polypeptide; Goodpasture antigen; OTTHUMP00000195044; TUMSTATIN; tumstatin
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)