

Product datasheet for **GA100900**

Human Collagen IV (COL4A4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	COL4A4 CRISPRa kit - CRISPR gene activation of human collagen type IV alpha 4 chain
Symbol:	COL4A4
UniProt ID:	P53420
Synonyms:	CA44
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)