


Product datasheet for **GA100913**

Human COL9A3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	COL9A3 CRISPRa kit - CRISPR gene activation of human collagen type IX alpha 3 chain
Symbol:	COL9A3
Components:	GA100913G1, COL9A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100913G2, COL9A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100913G3, COL9A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q14050
Synonyms:	DJ885L7.4.1; EDM3; IDD; MED
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)