

Product datasheet for **GA101022**

Human CSN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CSN3 CRISPRa kit - CRISPR gene activation of human casein kappa
Symbol:	CSN3
UniProt ID:	<u>P07498</u>
Synonyms:	CSN10; CSNK; KCA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)