

Product datasheet for **GA101079**

Human CX3CR1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CX3CR1 CRISPRa kit - CRISPR gene activation of human C-X3-C motif chemokine receptor 1
Symbol:	CX3CR1
UniProt ID:	P49238
Synonyms:	CCRL1; CMKBRL1; CMKDR1; GPR13; GPRV28; V28
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)