

Product datasheet for **GA101107**

Human Cytochrome P450 3A4 (CYP3A4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CYP3A4 CRISPRa kit - CRISPR gene activation of human cytochrome P450 family 3 subfamily A member 4
Symbol:	CYP3A4
UniProt ID:	P08684
Synonyms:	CP33; CP34; CYP3A; CYP3A3; CYP3A4; CYP3A5; HLP; NF-25; P450C3; P450PCN1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)