

Product datasheet for **GA101108**

Human Cytochrome P450 3A5 (CYP3A5) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CYP3A5 CRISPRa kit - CRISPR gene activation of human cytochrome P450 family 3 subfamily A member 5
Symbol:	CYP3A5
UniProt ID:	P20815
Synonyms:	CP35; CYP3A5; P450PCN3; PCN3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)