

## Product datasheet for **GA101133**

### Human DAGK (DGKQ) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	DGKQ CRISPRa kit - CRISPR gene activation of human diacylglycerol kinase theta
Symbol:	DGKQ
UniProt ID:	<a href="#">P52824</a>
Synonyms:	DAGK; DAGK4; DAGK7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)