

Product datasheet for **GA101137**

Human ZIP Kinase (DAPK3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DAPK3 CRISPRa kit - CRISPR gene activation of human death associated protein kinase 3
Symbol:	DAPK3
UniProt ID:	O43293
Synonyms:	DLK; ZIP; ZIPK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)