

Product datasheet for **GA101155**

Human **DCTN1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	DCTN1 CRISPRa kit - CRISPR gene activation of human dynactin subunit 1
Symbol:	DCTN1
UniProt ID:	Q14203
Synonyms:	DAP-150; DP-150; P135
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)