

Product datasheet for **GA101164**

Human DDOST activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | DDOST CRISPRa kit - CRISPR gene activation of human dolichyl-diphosphooligosaccharide--protein glycosyltransferase non-catalytic subunit |
| Symbol: | DDOST |
| UniProt ID: | P39656 |
| Synonyms: | AGER1; CDG1R; OKSWcl45; OST; OST48; WBP1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)