

Product datasheet for **GA101204**

Human DIO3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DIO3 CRISPRa kit - CRISPR gene activation of human iodothyronine deiodinase 3
Symbol:	DIO3
Components:	GA101204G1, DIO3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA101204G2, DIO3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA101204G3, DIO3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P55073</u>
Synonyms:	5DIII; D3; DIOIII; TXDI3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

