

## Product datasheet for **GA101205**

### Human **DKC1** activation kit by **CRISPRa**

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | DKC1 CRISPRa kit - CRISPR gene activation of human dyskerin pseudouridine synthase 1 |
| Symbol:       | DKC1   |
| UniProt ID:   | <a href="#">O60832</a>   |
| Synonyms:     | CBF5; DKC; DKCX; NAP57; NOLA4; XAP101  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)              |

#### Product images:



[View online >](#)