

Product datasheet for **GA101220**

Human Dystrophin (DMD) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | DMD CRISPRa kit - CRISPR gene activation of human dystrophin |
| Symbol: | DMD |
| UniProt ID: | P11532 |
| Synonyms: | BMD; CMD3B; DXS142; DXS164; DXS206; DXS230; DXS239; DXS268; DXS269; DXS270; DXS272; MRX85 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)