

Product datasheet for **GA101224**

Human Dystrophia myotonica protein kinase (DMPK) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | DMPK CRISPRa kit - CRISPR gene activation of human DM1 protein kinase |
| Symbol: | DMPK |
| UniProt ID: | Q09013 |
| Synonyms: | DM; DM1; DM1PK; DMK; MDPK; MT-PK |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:

