


Product datasheet for **GA101226**

Human dystrophin myotonia containing WD repeat motif (DMWD) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DMWD CRISPRa kit - CRISPR gene activation of human DM1 locus, WD repeat containing
Symbol:	DMWD
Components:	GA101226G1, dystrophin myotonia containing WD repeat motif gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA101226G2, dystrophin myotonia containing WD repeat motif gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA101226G3, dystrophin myotonia containing WD repeat motif gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q09019</u>
Synonyms:	D19S593E; DMR-N9; DMRN9; gene59
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)