

Product datasheet for **GA101228**

Human DNAH5 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | DNAH5 CRISPRa kit - CRISPR gene activation of human dynein axonemal heavy chain 5 |
| Symbol: | DNAH5 |
| UniProt ID: | Q8TE73 |
| Synonyms: | CILD3; DNAHC5; HL1; KTGNR; PCD |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)