

Product datasheet for **GA101238**

Human DYNC111 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DYNC111 CRISPRa kit - CRISPR gene activation of human dynein cytoplasmic 1 intermediate chain 1
Symbol:	DYNC111
Components:	GA101238G1, DYNC111 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA101238G2, DYNC111 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA101238G3, DYNC111 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14576</u>
Synonyms:	DNCI1; DNCIC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)