

Product datasheet for **GA101266**

Human Dopamine D2 Receptor (DRD2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DRD2 CRISPRa kit - CRISPR gene activation of human dopamine receptor D2
Symbol:	DRD2
UniProt ID:	P14416
Synonyms:	D2DR; D2R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)