

## Product datasheet for **GA101267**

### Human Dopamine Receptor D3 (DRD3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	DRD3 CRISPRa kit - CRISPR gene activation of human dopamine receptor D3
Symbol:	DRD3
UniProt ID:	<a href="#">P35462</a>
Synonyms:	D3DR; ETM1; FET1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)