

Product datasheet for **GA101268**

Human Dopamine Receptor D4 (DRD4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DRD4 CRISPRa kit - CRISPR gene activation of human dopamine receptor D4
Symbol:	DRD4
UniProt ID:	P21917
Synonyms:	D4DR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)