

Product datasheet for **GA101269**

Human Dopamine Receptor D5 (DRD5) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DRD5 CRISPRa kit - CRISPR gene activation of human dopamine receptor D5
Symbol:	DRD5
UniProt ID:	P21918
Synonyms:	DBDR; DRD1B; DRD1L2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)