

Product datasheet for **GA101271**

Human DRIL1 (ARID3A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARID3A CRISPRa kit - CRISPR gene activation of human AT-rich interaction domain 3A
Symbol:	ARID3A
Components:	GA101271G1, DRIL1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA101271G2, DRIL1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA101271G3, DRIL1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q99856
Synonyms:	BRIGHT; DRIL1; DRIL3; E2FBP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)